

## Daily Chamber Checklist (Sechrist)

Use a separate checklist for each chamber. Put a check on each line as completed unless otherwise indicated to enter a value or mark N/A if not applicable.

Week Of (MM/DD/YYYY) :

	MON	TUE	WED	THU	FRI
<b>Initials of person completing checklist</b>					
<b>Chamber Serial # _____</b>					
Oxygen supply pressure reading at alarm panel (50-90 psi) <b>Enter value</b> →					
Oxygen supply pressure reading at chamber (50-90 psi) <b>Enter value</b> →					
Air tank level (change below 500 psi) <b>Enter value</b> →					
Air tank regulator gauge set between 50-70 psi					
Air break equipment disinfected (demand valve only)					
Chamber covers completely removed and stored in a professional manner					
Chamber grounds are connected and without damage (back of chamber and patient grounding strap), check continuity with multimeter (less than 1 ohm Ω) <b>Enter value</b> →					
Chamber supply/vent hoses without kinks or damage					
Inspect acrylic chamber hull for scratches and/or crazing					
AC power connected and green LED lit					
Turn on entertainment					
With door open, rotate swing arm lever to closed position and turn master valve on. Check air flow and function of intercom (volume up and pick up handset)					
Inspect chamber door gasket for damage					
Inspect chamber controls for damage or loose knobs					
Inspect chamber interior and exterior for cleanliness					
<b>SHUTDOWN CHECKLIST</b>					
<b>Initials of person completing checklist</b>					
Turn off entertainment					
Turn the master valve to the OFF position					
Turn oxygen and air supply to chamber off (change air tank if below 500psi) <b>Enter value</b> →					
Clean chamber interior, exterior and cover					

Person Completing: \_\_\_\_\_ Initials \_\_\_\_\_

Person Completing: \_\_\_\_\_ Initials \_\_\_\_\_

Person Completing: \_\_\_\_\_ Initials \_\_\_\_\_

Person Completing: \_\_\_\_\_ Initials \_\_\_\_\_